

 Innovative Solutions. Proven Confidence.



AirHeel™

Specifically designed to treat plantar fasciitis, Achilles tendonitis, and heel pain. Two interconnected aircells located under the foot arch and in back of the heel, apply pulsating compression with every step, help reduce swelling and discomfort, and enhance circulation. The AirHeel is made of a breathable, lightweight fabric and is simple to apply and adjust with one strap.

DJO #	Sammons & Preston #	Shoe Size		Size
		Men	Women	
09AS	5629-52	≤ 7	≤ 8.5	S
09AM	9273-49	7-11.5	8.5-12.5	M
09AL	5629-53	11.5+	12.5+	L

Indications: Plantar fasciitis; Achilles tendonitis; Heel pain

Available from:



T 800.323.5547 F 800.547.4333
In Canada: T 800.665.9200 F 905.566.9207

In addition to the AirHeel, DJO offers a complete line of foot and ankle products.

ORDER ONLINE @ WWW.SAMMONSPRESTON.COM





PEER
REVIEWED
CLINICAL
SUPPORT

Plantar Fasciitis:

The Efficacy of a Pneumatic Compression Device in the Treatment of plantar Fasciitis

Kavros SJ, Mayo Clinic, Journal of Applied Biomechanics 404-413, 2005

- The Aircast AirHeel is efficacious in treating plantar fasciitis
- Patients with higher plantar fasciitis pain experience faster relief with the AirHeel than with a shoe insert

Achilles Tendonitis:

A Prospective Evaluation of the Pneumatic Achilles Wrap for treatment of Acute and Subacute Tendonitis

Mimms TT, Badekas A, Ottey DCK, Schon LC, Union Memorial Hospital, MD, 1999

- The Aircast AirHeel provides significant pain reduction for Achilles tendonitis

A Prospective Randomized Study Comparing the Therapeutic Effect of Eccentric Training, the AirHeel Brace, and a Combination of Both

Petersen W, Welp R, Rosenbaum D, American Journal of Sports Medicine, 2007

- The AirHeel is easier to use and as effective as eccentric training in the treatment of chronic Achilles tendinopathy
- After six week, pain levels decreased faster with the AirHeel than with the eccentric training

FEATURES &
BENEFITS

Features	Benefits
Intermittent pneumatic compression from two interconnected aircells under the foot arch and back of the heel	Swelling is minimized and circulation enhanced to provide rapid pain relief
Low profile sleeve	The AirHeel fits comfortably in most shoes
One-strap fitting	Application and adjustment of the AirHeel is easy even while wearing a shoe
Universal sizing	Small, Medium and Large sizes come in universal Left-Right, decreasing the inventory needed to carry

For more information, please contact your
Patterson Medical / Sammons Preston Sales Consultant



T 800.323.5547
F 800.547.4333

In Canada:
T 800.665.9200
F 905.566.9207

