

# Tumble Forms 2<sup>®</sup> Carrie<sup>®</sup> Seating System Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline 800-499-0285.

Account No. \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO#: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (    ) \_\_\_\_\_

## The flexible Carrie Seat can be used:

- At home — on a couch or in a chair. Tether strap secures seat in chair.
- On an airplane — elementary seat has been approved for use on an airplane.

### Carrie<sup>®</sup> Seat Features

- Shock absorbing foam
- Seamless covering, impervious to fluids and easy to clean
- Lightweight seat, easy to remove and position
- Height adjustable neck stabilizer and four-point adjustable harness

_____	Preschool Seat with Tray and Footrest	<b>4641S</b>	<b>\$1,836.20</b>
_____	Preschool Seat only	<b>9114</b>	<b>\$1,121.25</b>
_____	Preschool Footrest w/Bar	<b>9114S</b>	<b>\$459.95</b>
_____	Preschool Tray w/Bar	<b>9114T</b>	<b>\$409.95</b>
_____	Elementary Seat with Tray and Footrest	<b>4641M</b>	<b>\$1,974.95</b>
_____	Elementary Seat only	<b>9115</b>	<b>\$1,123.75</b>
_____	Elementary Footrest w/Bar	<b>9115S</b>	<b>\$481.20</b>
_____	Elementary Tray w/Bar	<b>9115T</b>	<b>\$431.20</b>
_____	Junior Seat with Tray and Footrest	<b>4641L</b>	<b>\$2,074.95</b>
_____	Junior Seat only	<b>9116</b>	<b>\$1,197.45</b>
_____	Junior Footrest w/Bar	<b>9116S</b>	<b>\$499.95</b>
_____	Junior Tray w/Bar	<b>9116T</b>	<b>\$448.70</b>
_____	Small Adult Seat with Tray and Footrest	<b>4641SA</b>	<b>\$2,437.45</b>
_____	Small Adult Seat only	<b>9117</b>	<b>\$1,343.70</b>
_____	Small Adult Footrest w/Bar	<b>9117S</b>	<b>\$771.20</b>
_____	Small Adult Tray w/Bar	<b>9117T</b>	<b>\$487.45</b>

**Tumble Forms 2<sup>®</sup>**  
 With Anti-Microbial Protection  
 Odor Resistance, Increased Durability



**Tumble Forms 2<sup>®</sup> Carrie<sup>®</sup> Seat Sizing Guide**

Dimensions	Preschool <b>4641S</b>	Elementary <b>4641M</b>	Junior <b>4641L</b>	Small Adult <b>4641SA</b>
Standing Height	30" - 38" (76 - 97 cm)	38" - 48" (97 - 122 cm)	48" - 58" (122 - 147 cm)	58" - 68" (147 - 173 cm)
Weight	20 - 40 lbs. (9 - 18 kg)	30 - 60 lbs. (14 - 27 kg)	50 - 100 lbs. (23 - 46 kg)	60 - 130 lbs. (27 - 59 kg)
Sitting Hip Width	6" - 8" (15 - 20 cm)	8" - 10" (20 - 25 cm)	10" - 13" (25 - 33 cm)	12" - 15" (31 - 41 cm)
Thigh Depth	8" - 11" (20 - 28 cm)	10" - 13" (25 - 33 cm)	12" - 15" (31 - 38 cm)	14" - 18" (36 - 46 cm)

Seat Dimensions	Pre-School <b>4641S</b>	Elementary <b>4641M</b>	Junior <b>4641L</b>	Small Adult <b>4641SA</b>
Inner Seat Height	22"	26"	30"	33"
Inner Seat Depth	8"	10"	13"	15½"
Inner Seat Width	8"	10"	12"	14"
Outside Seat Width	12"	14¼"	16¼"	18½"