

Leckey® Contoured Advance Seat Order Form

(page 1 of 2)



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

- Includes: Hi-Lo Base, push bar, cushioned armrests, individually adjustable footplates and multi-point pelvic strap.

Contoured Advance Seat

_____ Size 1	5668-78	\$4,106.20
_____ Size 2	5668-75	\$4,318.70
_____ Size 3, Manual Gas Lift	5668-76	\$4,499.95
_____ Size 3, Powered Lift	5668-77	\$5,184.95

Accessories (Size 1)

_____ Flexi Lateral Supports	5668-79	\$196.20
_____ Multi-Positional Contoured Headrest, Black Neoprene	5565-47	\$256.20
_____ Multi-Adjustable Contoured Headrest Upholstered	5649-32	\$244.95
_____ Flat Headrest	9244-01	\$269.95
_____ Laterals For Flat Headrest	9244-04	\$209.95
_____ Side Seat Pads	5649-34	\$223.70
_____ Activity Tray	5565-50	\$339.95
_____ Grab/Suspension Rail	5565-53	\$242.45
_____ Grab Posts	5565-55	\$181.20
_____ Padded Tray Insert (Grey)	5565-58	\$218.70
_____ Whitmyer Headrest Interface Bracket	5565-64	\$41.20
_____ Square Headrest Bracket (For headrests with a square post)	5535-79	\$53.70
_____ 1" Sandal Raiser (each)	9244-34	\$57.45
_____ Thoracic Strap	5565-66	\$142.45
_____ Chest Harness	5565-70	\$142.45
_____ Pommel	5668-99	\$167.45
_____ Rigid Laterals	0815-21-814	\$162.45

Accessories (Size 2)

_____ Flexi Lateral Supports	5656-92	\$207.45
_____ Multi-Positional Contoured Headrest, Black Neoprene	5565-48	\$268.70
_____ Multi-Adjustable Contoured Headrest Upholstered	5648-54	\$257.45
_____ Flat Headrest	9244-02	\$269.95
_____ Laterals For Flat Headrest	9244-05	\$213.70
_____ Side Seat Pads	5648-55	\$233.70
_____ Activity Tray	5565-51	\$359.95
_____ Grab/Suspension Rail	5565-54	\$242.45
_____ Grab Posts	5565-56	\$181.20
_____ Padded Tray Insert (Grey)	5565-59	\$218.70
_____ Whitmyer Headrest Interface Bracket	5564-64	\$41.20
_____ Square Headrest Bracket (For headrests with a square post)	5535-79	\$53.70
_____ 1" Sandal Raiser (each)	9244-34	\$57.45
_____ Thoracic Strap	5565-67	\$153.70
_____ Chest Harness	5565-71	\$153.70
_____ Pommel	5565-28	\$174.95
_____ Rigid Laterals	0815-21-822	\$162.45



Accessories (Size 3)

_____ Flexi Lateral Supports	5656-93	\$218.70
_____ Multi-Positional Contoured Headrest, Black Neoprene	5565-49	\$279.95
_____ Multi-Adjustable Contoured Headrest Upholstered	5648-60	\$271.20
_____ Flat Headrest	9244-03	\$279.95
_____ Laterals For Flat Headrest	9244-06	\$214.95
_____ Side Seat Pads	5648-61	\$246.20
_____ Activity Tray	5565-52	\$374.95
_____ Grab/Suspension Rail	5565-54	\$242.45
_____ Grab Posts	5565-56	\$181.20
_____ Padded Tray Insert (Grey)	5565-59	\$218.70
_____ Whitmyer Headrest Interface Bracket	5565-64	\$41.20
_____ Square Headrest Bracket (For headrests with a square post)	5535-79	\$53.70
_____ 1" Sandal Raiser (each)	9245-00	\$57.45
_____ Thoracic Strap	5565-68	\$164.95
_____ Chest Harness	5565-72	\$163.70
_____ Pommel	5668-98	\$184.95
_____ Rigid Laterals	0815-14-660	\$178.70

Leckey® Contoured Advance Seat Order Form

(page 2 of 2)




Contoured Advance Seat Patterns

Meadow



Size 1

Mondeo



Sizes 2 & 3

The Leckey Contoured Advance Seat comes with the fabric pattern indicated by the seat size and gray base.

Flipaway Kneecup Assembly

—— 2.625- 4" x 4" (65-100 x 100 mm), Size 1	5565-39	\$462.45
—— 2.625- 4" x 4.375" (65-100 x 110 mm), Size 2	5565-40	\$488.70
—— 2.625- 4" x 4.875" (65-100 x 120 mm), Size 3	5565-41	\$511.20
—— 2.625- 4" x 5.25" (65-100 x 130 mm), Size 4	5565-42	\$527.45

Sandals

—— 6.9" L x 3" W (175 x 75 mm), Size 1	5520-70	\$124.95
—— 8.3" L x 3.2" W (210 x 80 mm), Size 2	5520-71	\$134.95
—— 9.5" L x 3.8" W (240 x 95 mm), Size 3	5520-72	\$141.20

Age Range	Size 1 (1–5 years)	Size 2 (5–11 years)	Manual Size 3 (11–17 yrs)	Powered Size 3 (11–17 yrs)
Top of Seat to Floor	11"-26" (280-655 mm)	13"-26" (340-655 mm)	13"-26" (340-655 mm)	15"-27" (370-675mm)
Seat Depth	8"-12" (200-310 mm)	10"-15" (260-385 mm)	13"-19" (340-490 mm)	13"-19" (340-490 mm)
Seat Width	6"-12" (145-295 mm)	8"-14" (195-350 mm)	10"-16" (245-400 mm)	10"-16" (245-400 mm)
Seat to Sandal	7"-12" (190-300 mm)	10"-15" (245-390 mm)	13"-20" (320-520 mm)	13"-20" (320-520 mm)
Seat to Top Backrest	11"-14" (290-350 mm)	14"-16" (345-415 mm)	17"-22" (425-550 mm)	17"-22" (425-550 mm)
Back Width	6"-9" (140-240 mm)	7"-11" (180-290 mm)	9"-13" (220-340 mm)	9"-13" (220-340 mm)
Tray Height (from seat)	6"-10" (155-250 mm)	6"-10" (155-250 mm)	7"-13" (180-325 mm)	7"-13" (180-325 mm)
Max user weight	44 lbs (20 kg)	110 lbs (50 kg)	176 lbs (80 kg)	176 lbs (80 kg)
Tilt-in Space	Prone: 10°, Recline: 22.5°	Prone: 10°, Recline: 22.5°	Prone 10°, Recline: 22.5°	Prone 10°, Recline: 22.5°
Backrest Recline	Prone: 15°, Recline: 25°	Prone: 15°, Recline: 25°	Prone 15°, Recline: 25°	Prone 15°, Recline: 25°
Lower Back Support	-30° to 30°	-30° to 30°	-30° to +30°	-30° to +30°