

Leckey® Corner Sitter Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Features of the Corner Sitter

- Promotes upper body control
- Backrest is angle-adjustable from 20° recline to 15° prone
- Flat packed for easy storage
- Two built-in handles for easy portability

Corner Sitter

_____ Size 1	5649-35	\$904.95
_____ Size 2	5648-93	\$949.95
_____ Size 3	5648-95	\$999.95

Accessories

_____ Multi-Adjustable Contoured Head Support (Size 1)	5649-36	\$249.95
_____ Multi-Adjustable Contoured Head Support (Size 2)	5648-94	\$259.95
_____ Multi-Adjustable Contoured Head Support (Size 3)	5648-96	\$274.95
_____ Corner Sitter Tray Table (Fits all sizes)	9244-46	\$291.20



	Size 1	Size 2	Size 3
Age	1½-6	5-9	8-14
Overall User Height	45.9" (1165mm)	52.4" (1330mm)	64.2" (1630mm)
Seat Base Length	20.5" (520mm)	25.6" (650mm)	29" (740mm)
Seat Base Width	18.5" (470mm)	20.9" (530mm)	23.6" (600mm)
Height in Folded Position	7" (180mm)	7" (180mm)	7" (180mm)
Weight of Seat	33 lbs. (15kg)	35½ lbs. (16kg)	37 lbs. (17kg)
Backrest Height	11.8" (300mm)	15.7" (400mm)	19.7" (500mm)
Head Support	4.7" (120mm)	5.5" (140mm)	6.3" (160mm)
Head Support Adjustment	4" (100mm)	5.3" (135mm)	6.7" (170mm)
Pommel Adjustment	2.7" (70mm)	2.7" (70mm)	2.7" (70mm)
Max User Weight	66 lbs. (30kg)	88 lbs. (40kg)	110 lbs. (50kg)
Tray Height (from Ground)	8.8"-13.8" (225-350 mm)	8.8"-13.8" (225-350 mm)	8.8"-13.8" (225-350 mm)

Corner Sitter Patterns



The Corner Sitter comes with the fabric pattern indicated by the seat size.