

# Max II Positioning Chair Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at **800-547-4333** or call in your order to our Pediatric Helpline at **800-499-0285**.

Account No. \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO#: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (    ) \_\_\_\_\_

## Max II Tilt-n-Space Positioning Pushchairs Feature

- 45° tilt-in-space
- Full recline
- Supine position
- Suitable for children 1–12 years of age with significant growth adjustment
- Available in small or large size to accommodate a wide range of children
- Significant growth capacity for seat depth and seat width and back height
- Collapsible to compact size
- Meets crash tested ANSI/RESNA WC standards

Optional Sun/Rain Shield provides extra comfort and protection from the sun or rain.

## Max II Positioning Chair

\_\_\_\_ Max II Small    **5561-70**    **\$3,974.95**

## Max II Accessories

____ Sun/Rain Roof	<b>5561-72</b>	<b>\$349.95</b>
____ Sun/Rain Roof, Adjustable	<b>5561-73</b>	<b>\$624.95</b>
____ Headrest, Small	<b>5561-74</b>	<b>\$193.75</b>
____ Headrest, Large	<b>5561-75</b>	<b>\$206.25</b>
____ Abductor	<b>5561-76</b>	<b>\$64.95</b>
____ Vest Harness, Small	<b>5561-77</b>	<b>\$199.95</b>
____ Vest Harness, Large	<b>5561-78</b>	<b>\$199.95</b>
____ X-Harness, Small	<b>5561-79</b>	<b>\$228.75</b>
____ X-Harness, Large	<b>5561-80</b>	<b>\$228.75</b>
____ Pelvic Harness,	<b>5561-81</b>	<b>\$204.95</b>
____ Basket	<b>5561-88</b>	<b>\$118.75</b>



## Max II Specifications

Description	Seat Width	Seat Depth	Back Height	Compact Folding Length	Weight of Chair	User Wt. Capacity
Max II Small, Silver/gray	8" – 11¼"	8" – 12"	17" – 22"	31¾"	35 lbs.	77 lbs.