

# Special Tomato™ Multi-Positioning Seat Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston Rolyan at **800-547-4333** or call in your order to our Pediatric Helpline at **800-499-0285**.

Account No. \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO#: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (    ) \_\_\_\_\_

### Special Tomato MPS Seat Features:

- Individualized adjustment allows for many years of growth
- Seamless latex-free cushions are impermeable to fluids and soft, yet durable
- Cleans easily with common disinfectants.
- All Special Tomato MPS Seating Systems are FMVSS213 (Crash testing) Compliant

(special Tomato MPS Seats purchased before 10/03 are not suitable for use as a car seat.)



### Headrests. (Required) Fits both sizes of the MPS.

Aqua Blue	Lilac	Chocolate	Cherry	Lime	Description	
_____ 9288-62	_____ 9288-63	_____ 0814-51-814	_____ 0815-45-748	_____ 0815-45-904	Standard Headrest*	\$170.10
_____ 9288-64	_____ 9288-65	_____ 0814-51-822	_____ 0815-45-755	_____ 0815-45-912	Headrest with Laterals*	\$213.30

\*As of March 2006, a headrest must be purchased separately.

### Multi-Positioning Seat (MPS). Headrest not included. Includes Footrest.

Aqua Blue/Chocolate	Lilac/Chocolate	Chocolate/Chocolate	Cherry/Chocolate	Lime/Chocolate	Size of MPS Seat with Footrest	
_____ 9288-02	_____ 9288-24	_____ 0814-51-509	_____ 0815-45-524	_____ 0815-45-789	Small	\$1,785.55
_____ 9288-07	_____ 9288-29	_____ 0814-51-608	_____ 0815-45-631	_____ 0815-45-797	Small Extended	\$1,785.55
_____ 9288-13	_____ 9288-35	_____ 0814-51-574	_____ 0815-45-649	_____ 0815-45-805	Large	\$2,217.35
_____ 9288-18	_____ 9288-40	_____ 0814-51-665	_____ 0815-45-656	_____ 0815-45-813	Large Extended	\$2,217.35

### MPS Mobile Floor Kit. Includes MPS and Mobile Floor Base. Headrest not included. Includes Footrest.

Aqua Blue/Chocolate	Lilac/Chocolate	Chocolate/Chocolate	Cherry/Chocolate	Lime/Chocolate	Size of MPS Mobile Floor Kit	
_____ 9288-00	_____ 9288-22	_____ 0814-51-475	_____ 0815-45-664	_____ 0815-45-821	Small	\$2,678.35
_____ 9288-05	_____ 9288-27	_____ 0814-51-483	_____ 0815-45-672	_____ 0815-45-839	Small Extended	\$2,678.35
_____ 9288-11	_____ 9288-33	_____ 0814-51-616	_____ 0815-45-680	_____ 0815-45-847	Large	\$3,374.30
_____ 9288-16	_____ 9288-38	_____ 0814-51-640	_____ 0815-45-698	_____ 0815-45-854	Large Extended	\$3,374.30

### MPS Push Chair Kit. Includes MPS and Push Chair Base (9288-44). Headrest not included. Includes Footrest.

Aqua Blue/Chocolate	Lilac/Chocolate	Chocolate/Chocolate	Cherry/Chocolate	Lime/Chocolate	Size of MPS Push Chair Kit	
_____ 9288-01	_____ 9288-23	_____ 0814-51-491	_____ 0815-45-706	_____ 0815-45-862	Small	\$4,081.70
_____ 9288-06	_____ 9288-28	_____ 0814-51-533	_____ 0815-45-714	_____ 0815-45-870	Small Extended	\$4,081.70
_____ 9288-12	_____ 9288-34	_____ 0814-51-624	_____ 0815-45-722	_____ 0815-45-888	Large	\$4,579.55
_____ 9288-17	_____ 9288-39	_____ 0814-51-657	_____ 0815-45-730	_____ 0815-45-896	Large Extended	\$4,579.55

### MPS Car Seat. Includes MPS except for the footrest.

Aqua Blue/Chocolate	Cherry/Chocolate	Chocolate/Chocolate	Lilac/Chocolate	Lime/Chocolate		
_____ 9288-50	_____ 0815-46-027	_____ 0814-51-772	_____ 9288-54	_____ 0815-46-068	Small MPS Car Seat	\$1,099.75
_____ 9288-51	_____ 0815-46-035	_____ 0814-51-780	_____ 9288-55	_____ 0815-46-076	Small Extended MPS Car Seat	\$1,099.75
_____ 9288-52	_____ 0815-46-043	_____ 0814-51-798	_____ 9288-56	_____ 0815-46-084	Large MPS Car Seat	\$1,403.30
_____ 9288-53	_____ 0815-46-050	_____ 0814-51-806	_____ 9288-57	_____ 0815-46-092	Large Extended MPS Car Seat	\$1,403.30

# Special Tomato™ Multi-Positioning Seat Order Form (page 2 of 2)

five year  
warranty

Headrest sold  
separately



## Accessories

_____ MPS Push Chair Base	<b>9288-44</b>	<b>\$2,576.75</b>
_____ MPS Sun Shade for Push Chair Base	<b>9288-45</b>	<b>\$425.40</b>
_____ MPS Push Chair Travel Bag	<b>9288-47</b>	<b>\$425.40</b>
_____ Small Tomato Skin	<b>5507-46</b>	<b>\$120.60</b>
_____ Large Tomato Skin	<b>5507-47</b>	<b>\$128.30</b>
_____ MPS Pneumatic Tires Upgrade (Set of 4)	<b>0815-51-977</b>	<b>\$558.75</b>

## Mobile Floor Base

Built-in tray storage on back. Four swivel casters with rear locks. 10°-25° range of tilt-in space. 🗨️

Chocolate	Size	
_____ <b>9288-04</b>	Small	<b>\$924.50</b>
_____ <b>9288-15</b>	Large	<b>\$1,183.55</b>

## MPS Tray

Depth adjustable. 60° tilt adjustment. Tray easily attaches to back of Mobile Floor Base for storage.

Chocolate	Size	
_____ <b>9288-03</b>	Small	<b>\$412.70</b>
_____ <b>9288-14</b>	Large	<b>\$466.00</b>

## MPS Extended Seat Cushion

As your child outgrows the standard cushion, simply purchase the extended cushion.

Aqua Blue	Cherry	Chocolate	Lilac	Lime	Size	
_____ <b>9288-08</b>	_____ <b>0815-45-763</b>	_____ <b>0814-51-541</b>	_____ <b>9288-30</b>	_____ <b>0815-45-920</b>	Small	<b>\$283.15</b>
_____ <b>9288-19</b>	_____ <b>0815-45-771</b>	_____ <b>0814-51-558</b>	_____ <b>9288-41</b>	_____ <b>0815-45-938</b>	Large	<b>\$314.90</b>

## Footrests

Included with the MPS, MPS Mobile Floor Kit and MPS Push Chair Kit.

Chocolate	Description	
_____ <b>9288-58</b>	Small Footrest	<b>\$704.80</b>
_____ <b>9288-59</b>	Large Footrest	<b>\$835.60</b>



MPS Push Chair Base



Sun Shade



Tomato Skins



Mobile Floor Base



Extended Seat  
Cushion



MPS Tray

## Seating System Sizing Information

	Small		Large	
	Standard Seat	Extended Seat	Standard Seat	Extended Seat
Special Tomato Seats	Standard Seat	Extended Seat	Standard Seat	Extended Seat
Target User Age	2-5 years	5-8 years	8-11 years	11-14 years
Maximum User Wt.	80 lbs.	80 lbs.	150 lbs.	150 lbs.
Seat Weight	22 lbs.	24 lbs.	26 lbs.	28 lbs.
Seat & Floor Base Wt.	40 lbs.	42 lbs.	44 lbs.	46 lbs.
Seat Width	9"	10"	12"	13"
Seat Depth	8"	10"	13"	15"
Seat to Footrest	6" -16"	6" -16"	8" -18"	8" -18"
Top of Back Cushion	15" -19"	15" -19"	17" -21"	17" -21"
Top of Head Cushion	22" -27"	22" -27"	24" -29"	24" -29"
Hip Flex Angle	-15° to +15°	-15° to +15°	-15° to +15°	-15° to +15°
Floor to Seat Ht.	15"	15"	18"	18"
Tray	24"W x 17"D x 2"H	24"W x 17"D x 2"H	24"W x 17"D x 2"H	24"W x 17"D x 2"H
MPS Seat Only	15"W x 19"D x 42"H	15"W x 19"D x 42"H	19"W x 23"D x 47"H	19"W x 23"D x 47"H
MPS Floor Kit	19"W x 37"D x 37"H	19"W x 37"D x 37"H	22"W x 41"D x 42"H	22"W x 41"D x 42"H
MPS Push Chair Kit	25"W x 38"D x 42"H	25"W x 38"D x 42"H	25"W x 42"D x 47"H	25"W x 42"D x 47"H