



Leckey Squiggles Seating System Order Form

(Page 1 of 2)



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____
 Date: _____
 PO#: _____
 Buyer: _____
 Attention: _____

Ship to: _____
 Name: _____
 Address: _____
 City: _____
 Phone: () _____

Squiggles Seat Shell. Includes depth adjustable seat base, reclining backrest, lumbar positioning system, hip guides, and interface for Hi-Lo and Easy Base. Seat cushion/cover and pelvic harness are required

_____ 5619-70 **\$1,749.95**

Squiggles Seat Cushions/Covers. Includes backrest cushion, rear seat cushion, front seat cushion and sacral cushion. Requires pelvic harness.

_____ 5619-71 Green **\$374.95**
 _____ 5619-72 Orange **\$374.95**
 _____ 5619-73 Blue **\$374.95**
 _____ 5619-74 Pink **\$374.95**

Hi-Low Base. Manual gas assist with foot pedal activation.
 _____ 5619-75 **\$1,511.20**

Easy Base. With manual crank adjustment.
 _____ 5619-76 **\$1,371.20**

Push Bar Handle. Can be used with Hi-Low Base or Easy Base.
 _____ 5619-77 **\$218.70**

Convoid Safari Tilt Mobility Base. Includes Interface Bracket for Squiggles Seat.
 _____ 0814-33-564 **\$3,979.95**

Convoid EZ Rider Mobility Base. Includes Interface Bracket for Squiggles Seat.
 _____ 0814-33-648 **\$2,567.45**

Interface Bracket A (12"). (Required when attaching the seat to a wheelchair base except for the small Kimba, the Safari Tilt and EZ Rider bases.)
 _____ 5619-78 **\$258.70**

Trunk Harness Small, Black.
 _____ 5619-78-01 **\$223.70**

Trunk Harness Medium, Black.
 _____ 5619-78-02 **\$223.70**

Footrest Assembly.
 _____ 5624-56 **\$181.20**

Sandals Small. Includes sandal straps.*
 _____ 5520-70 **\$124.95**

Sandals Medium. Includes sandal straps.*
 _____ 5520-71 **\$134.95**

Sandals Raisers, 1". Single piece fits both sandals.
 _____ 5535-11 **\$68.70**

Tray.
 _____ 5619-79 **\$309.95**

Padded Tray Insert.
 _____ 5535-56 **\$249.95**

*Requires footrest assembly (5624-56)

Headrest Support Shell. (Requires Flat or Contoured Headrest Cushion.)
 _____ 5619-89 **\$172.50**

Square Headrest Bracket. For headrests with a square post.
 _____ 5535-78 **\$53.70**

Contoured Headrest Cushion. (Requires Headrest Support Shell 5619-89)
 _____ 5619-90 Green **\$87.45**
 _____ 5619-91 Orange **\$87.45**
 _____ 5619-92 Blue **\$87.45**
 _____ 5619-93 Pink **\$87.45**

Flat Headrest Cushion. (Requires Headrest Support Shell 5619-89)
 _____ 5619-94 Green **\$87.45**
 _____ 5619-95 Orange **\$87.45**
 _____ 5619-96 Blue **\$87.45**
 _____ 5619-97 Pink **\$87.45**

Flat Headrest Laterals. (For use with 5619-94 through 5619-97)
 _____ 5619-98 Green **\$259.95**
 _____ 5619-98-01 Orange **\$259.95**
 _____ 5619-98-02 Blue **\$259.95**
 _____ 5619-98-03 Pink **\$259.95**



Leckey Squiggles Seating System Order Form

(Page 2 of 2)

Rigid Lateral Trunk Supports. Includes covers. Pair.

_____ 5619-80 Green	\$256.20
_____ 5619-81 Orange	\$256.20
_____ 5619-82 Blue	\$256.20
_____ 5619-83 Pink	\$256.20

Flip Away Lateral with Cover.

_____ 0815-14-728 Green	\$437.45
_____ 0815-14-736 Orange	\$437.45
_____ 0815-14-751 Pink	\$437.45
_____ 0815-14-744 Blue	\$437.45

Chest Support Harness. (Requires Rigid Lateral Trunk Supports)

_____ 5619-84 Green	\$144.95
_____ 5619-85 Orange	\$144.95
_____ 5619-86 Blue	\$144.95
_____ 5619-87 Pink	\$144.95

Pelvic Harness Spacer Pads.

_____ 5619-99	\$69.95
---------------	---------

(For use with 5619-99-01 through 5619-99-04. Provides extra cushioning to pelvic harness.)

Small Pelvic Harness.

_____ 5619-99-01 Green	\$124.95
_____ 5619-99-02 Orange	\$124.95
_____ 5619-99-03 Blue	\$124.95
_____ 5619-99-04 Pink	\$124.95

Medium Pelvic Harness.

_____ 0815-22-077 Green	\$122.45
_____ 0815-22-085 Orange	\$122.45
_____ 0815-22-093 Blue	\$122.45
_____ 0815-22-101 Pink	\$122.45

Seat Pommel.

_____ 5619-88	\$87.45
---------------	---------

1" Back Pad and Rigid Laterals.

Reduces seat depth by 1". Includes special laterals for the thicker back pad.

_____ 5535-46 Green	\$484.95
_____ 5535-47 Orange	\$484.95
_____ 5535-48 Blue	\$484.95
_____ 5535-49 Pink	\$484.95

Grab Rail.

_____ 5520-79	\$226.20
---------------	----------

Squiggles color swatches



Squiggles Seat sold separately.



Shown on Convoid Safari Tilt Mobility Base.



Squiggles Seat Specifications:

Age (approx):	1-5
Max User Weight:	48.4 lbs. (22kg)
User Height:	Min: 29.5" (750mm) Max: 43.3" (1100mm)
Chest Width between laterals:	Min: 6.25" (160mm) Max: 8.75" (220mm)
Seat Width between hip guides:	Min: 6.25" (160mm) Max: 8.75" (220mm)
Seat Depth:	Min: 7.5" (190mm) Max: 11.4" (290mm)
Seat Depth with 1" Back Pad:	Min: 6.5" (216mm) Max: 10.4" (315mm)
Knee Width between pommel and side of seat:	3.5" (90mm)
Backrest Height:	Min: 12.6" (320mm) Max: 15" (380mm)
Backrest Angle:	Prone 10°; Recline 25°
Seat to Footplate: (soul of foot to back of knee)	Min: 6.3" (160mm) Max: 10.6" (270mm)
Footplate Angle Adjustment:	10° Plantarflexion/Dorsiflexion
Top of Seat to Floor (Hi Low Base):	Min: 10.6" (270mm) Max: 25.6" (650mm)
Top of Seat to Floor (Easy Base):	Min: 11" (280mm) Max: 15" (380mm)
Tray Height:	Min: 4.5" (115mm) Max: 6.75" (170mm)
Tray Size:	18.1" x 22" (460 x 560mm)
Sandal Size:	Small 6.9"L x 3"W (175-75mm) Medium 8.3"L x 3.2"W (210-80mm)
Seating Unit Weight:	13.2 lbs. (6kg)
Minimum dimensions for storage:	14.6"L x 14.6"W x 15.7"H (370 x 370 x 400mm)

Hi Low Base (excluding Push Bar)

Weight:	19.8 lbs. (9kg)
Min. dimensions:	25.2"L x 21.3"W x 12.4"H (640 x 540 x 315mm)

Easy Base (excluding Push Bar)

Weight:	9.9 lbs. (4.5kg)
Min. dimensions:	24.8"L x 20.7"W x 10.4"H (630 x 525 x 265mm)

Convoid Safari Tilt Mobility Base:

Weight:	31.4 lbs. (14.2kg)
Overall Width:	22.5" (572mm)
Overall Depth:	31" (787mm)
Handle Height:	44.5" (1130mm)
Wheel Diameter Front:	6" (152mm)
Wheel Diameter Rear:	11" (279mm)
Tilt Range:	5-45 degrees
Weight Capacity:	48 lbs (21.8kg)
Folded:	50"L x 14"W x 14"D (1270 x 360 x 360mm)

Convoid EZ Rider Mobility Base Specifications

Weight:	20.4 lbs. (9.3kg)
Overall Width:	18.88" (480mm)
Overall Depth:	29.25" (743mm)
Handle Height:	36.5" (927mm)
Wheel Diameter, Front:	6" (152mm)
Wheel Diameter, Rear:	8" (203mm)
Recline Angle:	10°
Weight Capacity:	48 lbs. (21.8kg)
Folded:	43"L x 13.5"W x 14.5"D (1090 x 340 x 370mm)